Cranibas

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The ARTZT vitality Cranibas is based on a therapeutic concept, thanks to which the area surrounding the articulation between the cervical vertebrae and the skull can be gently stretched and relaxed, simply by correctly resting on it. The intrinsic weight of the head causes gentle stretching and relaxation of the connective tissue structures in the area of the articulation between the cervical vertebrae and the skull (muscles, ligaments and fascia). Furthermore, metabolic processes and elasticity can be influenced and the mobility of the articulation between the cervical vertebrae and the skull can be returned to normal.

General instructions

- » Only carry out these exercises in consultation with a medical professional.
- » The ARTZT vitality Cranibas should be kept away from the eyes to avoid possible injuries by the cusps.
- » The ARTZT vitality Cranibas is not a toy and should be kept out of the reach of children.

Please note: The ARTZT vitality Cranibas therapeutic pad is a therapeutic aid and is not a substitute for a doctor's appointment, nor for osteopathic or physiotherapeutic consulting/treatment. Contact your doctor or physiotherapist before using the ARTZT vitality Cranibas.

The ARTZT vitality Cranibas is a product derived from practical experience, not from scientific studies. The following descriptions are based solely on experience with patients in dialogue with physiotherapists and osteopaths and do not claim to have general validity.

Since the treatment result of the Cranibas depends to a large extent on the regularity of its use, therapy should not be interrupted even during pain-free periods. If its use is painful at first, the recommended initial period of use of five minutes should not be exceeded. However, a feeling of pressure and tension and possibly slight pain is quite normal and an expression of a specific stimulus to the region below the back of the head, which is heavily stressed in everyday life.

Function

The various angles and elevations of the Cranibas therapy pad serve to support the head and the area between the cervical vertebrae and the base of the skull. Due to the special arrangement of the cusps, slight pressure and pull is exerted in the area below the back of the head. This mechanical stimulus can affect tension in the connective tissue structures in this region (muscles, ligaments and fascia).

Material and characteristics

The ARTZT vitality Cranibas is made of PU foam with a skin-friendly surface. All materials used are odourless and well tolerated by the skin.

Available in two different versions:

Home: Velvety flock for private use.

Professional: Smooth surface for use in therapy and practice.

Cleaning and maintenance

The surface of the ARTZT vitality Cranibas Home and Professional can be optimally treated with PU-compatible surface disinfectant or cleaned with a damp chamois and some household soap.

Put on

Sit on a stable, even surface (e.g. floor or gymnastic mat). Place both cusps of the ARTZT vitality Cranibas on the lower edge of the back of the head, at the level of the base of the skull, or approximately at the hairline. The large hollow should be facing upwards.

Lie down

Place the large hollow of the ARTZT vitality Cranibas on the back of your head. Tilt your head backwards until the ARTZT vitality Cranibas touches your neck. Maintaining this position, settle carefully into a supine position. If possible, stretch out your legs. The ARTZT vitality Cranibas must rest completely on the undersurface or on the floor.

Relax

Let ARTZT vitality Cranibas work. A sensation of gentle pressure build-up is entirely normal and a sign of gentle stretching at the base of the skull.

Stay in this position for a maximum of 15 minutes!

Usage period

It should initially be used for periods not exceeding five minutes and then treatment time gradually increased. It should not be used for periods of more than 15 minutes. Shorter and more frequent use (e.g. morning and evening for 10 minutes each) is preferable to longer periods of use. In the event of any discomfort during use, this must be discontinued and a medical professional must be consulted.

Intensifying the effect

- » Move the eyes, e.g. in circular movements.
- » Gentle, slow circular head movements mobilise the articulation between cervical vertebrae and the skull.
- » Slight lifting and lowering of the feet with outstretched legs intensifies the stretching.
- » Heat application prior to use (e.g. by placing a hot-water bottle on the neck) results in faster relaxation of the corresponding structures.
- » For optimal mobilisation of the articulations between the head and cervical vertebrae, slightly stretch your head in the direction of rotation.

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©2018 Ludwig Artzt GmbH. All rights reserved. Printed in Germany.
Ludwig Artzt GmbH
©2018 Ludwig Artzt GmbH. Alle Rechte geschützt. Gedruckt in Deutschland.